


DNA Stain in KRISHAN BUFFER (PI solution 0.05 mg/ml)

From: SOP RPCI Department of Flow and Image
Cytometry – 2013



Reagents:

- 0.1% Sodium Citrate
- 0.2 mg/ml RNase (Worthington # 5356965).
Procedure note: Because RNase is difficult to weigh out, a stock solution of 2.5 mg/ml is made up and stored at 4°C.
- 0.05 mg/ml Propidium Iodide (Sigma #P-4170)
- 0.2% NP 40 (or Igepal, Sigma #I-3021)
- 1N HCl
- Deionized H₂O.

Procedure:

1. Place 200 ml of deionized H₂O into dark bottle
2. Add:
 - 0.2 grams of Sodium Citrate
 - 1.6 ml 2.5 mg/ml RNase solution
 - 10 mg Propidium Iodide
 - 400 µl of NP 40
 - 1 drop 1N HCl(all volumes may be divided in half to make up only 100 ml of stain)
3. Store at 4 °C (shelf life 3-4 weeks).